Heating and cooling thermostats with temperatures from -90 up to 300 °C for professional use in research, application engineering and production









Application examples

- Temperature control for chemical synthesis
- Tests on electronic components at different temperatures
- Temperature control of measuring structures in process technology
- Heating and cooling of glass reactors

Proline Kryomats

Application examples

- Constant temperatures in the notch bending test and drop test
- Changing temperatures when determining the pour point, Brookfield test of lubricants and test of slide bearings



Intuitive operation, ultra high and low temperatures

LAUDA Proline thermostats are our solution for high performance and reliable temperature regulation. With their broad temperature range they fulfill high requirements. LAUDA Proline thermostats are available in two versions: in the basic version with the Master control head, and as a Command version with a removable control unit for enhanced ease of operation. Master version devices can be retrofitted with the Command remote control, which are simply connected to the control head. The thermostat automatically recognises and controls all newly installed moduls.

Your advantages at a glance

+	The Proline advantages	Your benefits
	 Master or Command version 52 different devices Simple retrofitting from Master to Command version 	 The right solution for every application Subsequent extension or adaptation to changing application requirements
	Graphical user guidanceAdaptive control on cooling thermostats	 Easy and intuitive operation Saves time-consuming calculation of control parameters
	 Patented SmartCool system PowerAdapt system for adjustment of the power consumption 	 Up to 75 percent energy saving with digital cooling management Use of the maximum available output from the power supply system
	 Two insert ports can be combined with five different interface modules Easy distribution of the pump flow by means of bypass valve Pump connections on the side and rear 	 High level of flexibility for the user allowing for broad range of system integration Simultaneous connection of two external applications Flexible connection of external appli- cations from different sides
	 High-performance pressure-suction pump (Varioflex pump) with eight pump levels Up to 3.5 kW (230 V) heating power – even on all cooling thermostats via SmartCool system 	 Suitable for internal and external applications Adaptation of the pump power to the respective application Rapid heating achieved

Proline Master control head

The Proline Master devices are designed with high thermostating accuracy and reliability for all applications from -90 up to 300 °C where operating parameters are not changed or modified frequently. They have all the basic features and safety functions required for professional thermostating during continuous use. A modular structure and bus technology have created an instrument capable of extending its function and performance as the application requires.





- Easy-to-read green LED display
- Convenient setting of set-temperature and Varioflex pump via three operation buttons
- Indicator lights for heating, cooling, external control and alarm
- Resolution of indication 0.01 °C, setting resolution selectable 0.1 or 0.01 °C
- Selectable operating temperature range and additional button for overtemperature protection setting
- External temperature control via Pt100
- Optical and audible alarm function
- Simple temperature probe calibration
- Integrated mains network safety device
- Start mode control (automatic or manual)
- Two slots for interface modules
- LAUDA Wintherm Plus control software via RS 232/485 interface (optional)



Easy replacement of interface modules



Alarm message for malfunction



Upgradable to Command version

Proline Command control head

The Command control heads are the top models of the LAUDA Proline. The highly-efficient programer fulfills all the requirements of complex thermostating processes – with real-time function. It offers the utmost in user-friendliness and optimum functionality, e.g. for an industrial testing lab. The simple menu-driven operation and the easy editing of test programs allow for quickly changing thermostating tasks. The Command remote control is removable and can easily be used with cable up to 50 m. Comprehensive basic equipment as with the Proline Master range.





Basic equipment as Proline Master, plus:


- 10-key console for setpoint adjustment of temperature
- Programer with real-time clock, 150 temperature/time segments, for use in up to 5 programs, editable segments with loop and tolerance band function
- High resolution, back-lit, graphic LCD display with various display possibilities
- Detachable Command remote control for use with cable up to 50 m
- Eight freely selectable fixed temperatures with memory function
- Resolution of actual value display up to 0.001 °C
- RS 232/485 interface for LAUDA Wintherm Plus software
- Menu guidance in German, English, French and Spanish



An opto-decoupled RS 232/485 interface is integrated as standard

Pump	Level
Settings	Calibration
Graph	Default Settings
Clock	Resolution
Programmer	Device Status
Interfaces	Keyboard
Control	Basic Settings
Limits	
Pump Menu	End Tset Tfix

Drop-down menus make settings easy. Available in four languages.




Proline Heating thermostats with Master control head up to 19 liters

The heating thermostats of the LAUDA Proline with Master control head do not only convince because of their compact construction. The high heater power of 3.5 kW (230 V), two interfaces for various modules, a cooling coil fitted as a standard feature, and an integrated external control – these features make them particularly useful for users who require flexible thermostating operations while only rarely needing to adjust the settings.



Heating thermostat P 18



Temperature range 30...300 °C

Included accessories

Bath cover \cdot 2 nipples and 4 closing plugs for pump connections \cdot 2 nipples for cooling coil

Additional accessories

454 mm

Constant level device (for P 8) · automatic filling device · through-flow cooler · reverse flow protection · tubes · solenoid valve for cooling water control · high-temperature cooler (water) · Interface modules: analog, RS 232/485, contact, Profibus module

574 mm



All technical data on page 92 and following Other power supply variants on page 103

Other power suppl	ly variants on page 103				
Technical features		P 5	P 8	P 12	P 18
Working temperature range*	°C	35300	35300	30300	30300
Temperature stability	±Κ	0.01	0.01	0.01	0.01
Heater power	kW	3.5	3.5	3.5	3.5
Pump pressure max.	bar	0.7	0.7	1.1**	0.7
Pump suction max.	bar	0.4	0.4	-	0.4
Pump flow (pressure) max.	L/min	25	25	32**	25
Pump flow (suction) max.	L/min	23	23	-	23
Bath volume	L	3.55.5	5.58	6.513.5	12.519
Bath opening/Bath depth	mm	150x50/200	150x150/200	150x150/320	300x200/200
Cat. No. 230 V; 50/60 Hz		LCB 0708	LCB 0710	LCB 0716***	LCB 0712

Applications Advantages Control heads Devices Accessories

* Working temperature range with water cooling 20...300 °C ** Pressure pump only, pump characteristics see page 42 *** Instead of pressure and suction pump equipped with increased output


Proline Heating thermostats with Master control head up to 53 liters

The LAUDA Proline P 26, P 40 and P 50 heating thermostats are distinguished by particularly large-volume baths. All the below models are equipped with a Varioflex pump and cover the temperature range from 30 up to 300 °C. These stainless steel baths are ideally suited to direct thermostating inside the bath. The P 40 is particularly suitable for thermostating applications needing a large submersion depth. The P 26 and P 50 models with their wide baths, allow long or bulky test pieces to be placed in the bath or even enable a number of test pieces to be positioned alongside each other, for simultaneous testing.

A circulation chamber on the P 40 and P 50 ensures good mixing in the bath and thus guarantees good temperature homogeneity, despite the large bath vessel.



Heating curves Heat transfer liquid: Ultra 300, bath closed



Temperature range 30...300 °C

Included accessories

Bath cover (only P 26) \cdot 2 nipples and 4 closing plugs for pump connections \cdot 2 nipples for cooling coil


Additional accessories

Bath cover · automatic filling device · through-flow cooler · reverse flow protection · tubes · solenoid valve for cooling water control · high-temperature cooler (water) · rising platform · Interface modules: analog, RS 232/485, contact, Profibus module



All technical data on page 92 and following Other power supply variants on page 103





* Working temperature range with water cooling 20...300 °C ** Max. temperature only achieved with closed bath cover


454 mm

Proline Heating thermostats with Command control head up to 19 liters

The Proline heating thermostats with Command control head (C) impress through an expanded scope of functions. Alongside a graphic LCD display, which enables current values to be displayed up to 0.001 °C resolution, an easily editable and convenient programmer with storage possibilities is available. The standard RS 232/485 interface enables communication with a computer. Work flexibly with Command: The Command remote control can be quickly and easily detached from the thermostat.



Heating thermostat P 18 C



Temperature range 30...300 °C

20

Included accessories

Bath cover \cdot 2 nipples and 4 closing plugs for pump connections \cdot 2 nipples for cooling coil

40

60

min Heating time

Additional accessories

Constant level device (for P 8 C) · automatic filling device · through-flow cooler · reverse flow protection · tubes · solenoid valve for cooling water control · hightemperature cooler (water) · Interface modules: analog, RS 232/485, contact, Profibus module

All technical data on page 92 and following Other power supply variants on page 103

Other power supp	iy vananis on page 105				
Technical features		P 5 C	P 8 C	P 12 C	P 18 C
Working temperature range*	°C	35300	35300	30300	30300
Temperature stability	±Κ	0.01	0.01	0.01	0.01
Heater power	kW	3.5	3.5	3.5	3.5
Pump pressure max.	bar	0.7	0.7	1.1**	0.7
Pump suction max.	bar	0.4	0.4	-	0.4
Pump flow (pressure) max.	L/min	25	25	32**	25
Pump flow (suction) max.	L/min	23	23	-	23
Bath volume	L	3.55.5	5.58	6.513.5	12.519
Bath opening/Bath depth	mm	150x50/200	150x150/200	150x150/320	300x200/200
Cat. No. 230 V; 50/60 Hz		LCB 0709	LCB 0711	LCB 0717***	LCB 0713

* Working temperature range with water cooling 20...300 °C ** Pressure pump only, pump characteristics see page 42 *** Instead of pressure and suction pump equipped with increased output


Proline Heating thermostats with Command control head up to 53 liters

In order to enhance ease of use even further, P 26 C, P 40 C and P 50 C heating thermostats with large baths are also available with the Command control head which allows complex thermostating functions, particularly those with internal thermostating processes, to be easily mastered with the aid of an intuitive operation guidance system and the ability to edit programs rapidly.

A circulation chamber on the P 40 and P 50 ensures good mixing in the bath and thus guarantees good temperature homogeneity, despite the large bath vessel.



Heating thermostat P 40 C



Temperature range 30...300 °C

Included accessories

Bath cover (only P 26 C) \cdot 2 nipples and 4 closing plugs for pump connections \cdot 2 nipples for cooling coil

Additional accessories

Automatic filling device · bath cover · through-flow cooler · reverse flow protection · tubes · solenoid valve for cooling water control · high-temperature cooler (water) · rising platform · Interface modules: analog, RS 232/485, contact, Profibus module



All technical data on page 92 and following Other power supply variants on page 103





Technical features		P 26 C	P 40 C	P 50 C
Working temperature range*	°C	30300	30300**	30300**
Temperature stability	±Κ	0.01	0.01	0.01
Heater power	kW	3.5	3.5	3.5
Pump pressure max.	bar	0.7	0.7	0.7
Pump suction max.	bar	0.4	0.4	0.4
Pump flow (pressure) max.	L/min	25	25	25
Pump flow (suction) max.	L/min	23	23	23
Bath volume	L	1827	3037	3553
Bath opening/Bath depth	mm	300x350/200	250x270/450	300x750/200
Cat. No. 230 V; 50/60 Hz		LCB 0715	LCB 0729	LCB 0731

* Working temperature range with water cooling 20...300 °C ** Max. temperature achieved only with closed bath cover

Proline Viscothermostats

LAUDA viscothermostats are optimized for directly observing inserted objects. The temporal and spatial temperature stability required for precisely determining the viscosity is guaranteed for the full temperature range. As such, they are ideal for use with the fully automated LAUDA PVS or iVisc viscometers. Thanks to the double-chamber principle, a constant liquid level in the measuring room is guaranteed regardless of the rate and temperature. The PVL models are equipped with five layers of insulating glass and by connecting a DLK 45 through-flow cooler or Proline RP 890 cooling thermostat are suited to low-temperature measurements down to -40 or -60 °C.




Viscothermostat PV 24 C




All technical data on page 92 and following Other power supply variants on page 103

Pump characteristics Heat transfer liquid: Water



Heating curves Heat transfer liquid: Therm 240, bath closed



Temperature range 30...230 °C

Included accessories

2 nipples and 4 closing plugs for pump connections · 2 nipples for cooling coil

Additional accessories Window heating system – PVL 15 (C), PVL 24 (C) only \cdot solenoid valve for cooling water \cdot additional cooler \cdot Command

remote control

Technical features		PV 15/PV 15 C	PV 24/PV 24 C	PV 36/PV 36 C	PVL 15/PVL 15 C	PVL 24/PVL 24 C
Working temperature range	°C	30230	30230	30230	30100	30100
Temperature stability	±Κ	0.01	0.01	0.01	0.01	0.01
Heater power	kW	3.5	3.5	3.5	3.5	3.5
Pump pressure max.	bar	0.8	0.8	0.8	0.8	0.8
Pump suction max.	bar	-	-	-	-	-
Pump flow (pressure) max.	L/min	25	25	25	25	25
Pump flow (suction) max.	L/min	-	-	-	-	-
Bath volume	L	1115	1924	2836	1115	1924
Bath opening/Bath depth	mm	230x135/320	405x135/320	585x135/320	230x135/320	405x135/320
Glass pane size	mm	149x230	326x230	506x230	149x230	326x230
Cat. No. Master 230 V; 50/60 Hz		LCD 0276	LCD 0278	LCD 0280	LCD 0282	LCD 0284
Cat. No. Command 230 V; 50/60 Hz		LCD 0277	LCD 0279	LCD 0281	LCD 0283	LCD 0285

Proline Bridge thermostats

LAUDA Proline bridge thermostats are available in two versions with different pump models and immersion depths. The PB models have a pressure/suction pump and require a bath depth of 200 mm, while the PBD models have a more powerful pressure pump (D) and require a bath with a depth of 320 mm. In addition, both series of models differ in the selected control head: Master or Command (C). Through variably extendable telescopic rods, all models can be attached without problem to baths with a width from 310 mm up to 550 mm.





Bridge thermostat PBD C

- Bath not included in scope of delivery -




All technical data on page 92 and following Other power supply variants on page 103





Pump characteristics for PBD and PBD C P 12 and P 12 C, Heat transfer liquid: Water



Temperature range 30...300 °C

Included accessories

2 nipples and 4 closing plugs for pump connections \cdot telescopic rods

Additional accessories

Automatic filling device · water bath


Interface modules: analog, RS 232/485, contact, Profibus module

Technical features		PB/PB C	PBD/PBD C	
Working temperature range	°C	30300	30300	
Operating temperature range	°C	-30*300	-30*300	
Temperature stability	±Κ	0.01	0.01	
Heater power	kW	3.5	3.5	
Pump pressure max.	bar	0.7	1.1	
Pump suction max.	bar	0.4	_	
Pump flow (pressure) max.	L/min	25	32	
Pump flow (suction) max.	L/min	23	_	
Bath volume up to approx.	L	80	80	
Bath opening	mm	Telescopic rods can be exte	ended for bath widths 310550	
Bath depth min.	mm	200	320	
Cat. No. Master 230 V; 50/60 Hz		LCG 0090	LCG 0092	
Cat. No. Command 230 V; 50/60 H	łz	LCG 0091	LCG 0093	
* Only achievable with LAUDA through-flow co	oler			

43

Proline Cooling thermostats with Master control head up to 8 liters

The Proline RP 845, RP 855, RP 870 and RP 890 cooling thermostats stand out above all for their compact dimensions and small footprint. With a cooling capacity of 1.6 kW at 20 °C, the RP 855 has a particularly highperformance design. The RP 890 low-temperature device enables you to reach temperatures down to -90 °C. The standard integrated electricalheated bath cover set prevents icing as a result of condensation and humidity on all Proline cooling thermostats with a temperature range down to -90 °C.



LCK 1885

LCK 1893

LCK 1895

LCK 1897

*

Bath temperature

°C

20

0

Cooling curves Heat transfer liquid: Ethanol, bath closed

① RP 845


@ RP 3530 3 RP 1840

* Working temperature range is equal to ACC range

Cat. No. 230 V; 50 Hz

Proline Cooling thermostats with Master control head up to 35 liters

The Proline RP 1290, RP 1840, RP 1845 and RP 3530 cooling thermostats differ in terms of bath volume, achievable cooling capacity, and working temperature range. With bath capacities up to 35 liters, the RP 3530 provides a particularly large bath volume and the RP 1845 provides a particularly high cooling capacity of 1.6 kW. To prevent icing the RP 1290 is equipped with an electrical-heated bath cover set.



Temperature range -88...200 °C

Included accessories Bath cover · 2 nipples and 4 closing plugs for pump connections · electrical-heated bath cover set for RP 129 Additional accessories Reverse flow protection · automatic filling device · tubes Interface modules: analog, RS 232/485, contact, Profibus module · set of castors (RP 1290 to RP 1845)				illing device · tubes · 85, contact, Profibus	
Cooling thermostat RP 1845		Т			T
All technical data on and following Other power supply va		04	688 mm	770 mm	740 mm
Technical features		RP 1290	RP 1840	RP 1845	RP 3530
Working temperature range*	°C	-88200	-40200	-50200	-35200
Temperature stability	±Κ	0.02	0.01	0.01	0.02
Heater power	kW	3.5	3.5	3.5	3.5
Cooling output at 20 °C	kW	1.1	0.9	1.6	0.9
Pump pressure max.	bar	0.7	0.7	0.7	0.7
Pump suction max.	bar	0.4	0.4	0.4	0.4
Pump flow (pressure) max.	L/min	25	25	25	25
Pump flow (suction) max	L/min	23	23	23	23
Bath volume	L	813.5	12.519	12.519	2335
Bath opening/depth	mm	300x150/200	300x200/200	300x200/200	300x350/250
Cat. No. 230 V; 50 Hz		LCK 1899	LCK 1887	LCK 1891	LCK 1889

₩

* Working temperature range is equal to ACC range

Proline Cooling thermostats with Command control head up to 8 liters

The SmartCool system – an energy-saving, digital cooling management system – ensures that every temperature is run with the correct cooling capacity. It increases or reduces the cooling according to application requirements. The advantages are particularly effective for programmer operation and temperature ramping.

The Proline cooling thermostats with the Command control head (C) have a convincing extended range of functions. At 20 $^{\circ}$ C, RP 855 C has a particularly high cooling capacity of 1.6 kW. RP 890 C and RP 1290 C are designed for especially low temperatures. They differ in terms of bath volume and have an electrical-heated bath cover set as a standard feature.



L/min

L


mm

23

5.5...8

150x150/200

LCK 1886



Temperature range -90...200 °C


*

Included accessories

Bath cover \cdot 2 nipples and 4 closing plugs for pump connections \cdot electrical-heated bath cover set for RP 890 C

Additional accessories

Constant level device for RP 845 C und RP 855 C \cdot reverse flow protection \cdot automatic filling device \cdot tubes \cdot electricalheated bath cover set for RP 855 C and RP 870 C (only ex works) \cdot Interface modules: analog, RS 232/485, contact, Profibus module \cdot set of castors (RP 855 C to RP 890 C)



23

5.5...8

150x150/200

LCK 1894

23

5.5...8

150x150/200

LCK 1896

23

5.5...8

150x150/200

LCK 1898

Applications Advantages Control heads Devices Accessories

Cat. No. 230 V; 50 Hz

Pump flow (suction) max.

Bath volume

Bath opening/depth

* Working temperature range is equal to ACC range

Heating curves Heat transfer liquid: Ultra 300, bath closed

4

① RP 855 C

RP 845 C

RP 870 C

RP 890 C

② RP 1290 C

③ RP 1840 C

④ RP 3530 C

RP 1845 C

Proline Cooling thermostats with Command control head up to 35 liters

Thanks to their various capacity ranges and filling volumes, the Proline cooling thermostats which make up the Command range skilfully adapt to your requirements. The RP 1845 C works at a temperature range between -50 and 200 °C and, at 20 °C, has a cooling capacity of 1.6 kW. The RP 3530 C has a particularly large bath for internal sample thermostating. The RP 1290 C comes with an integrated electrical-heated bath cover set as standard.



*

Bath temperature

23

1

°C

200

150

100

50

* Working temperature range is equal to ACC range

LAUDA Proline Kryomats

Extra powerful cooling thermostats for bath applications from -90 up to 200 °C LAUDA Proline Kryomats



The **Proline Kryomats** are floor-standing, low temperature thermostats suitable for a wide variety of applications. They never fail to impress through their compact design and high cooling capacities, especially at low temperatures. All Proline Kryomats are fitted with the Command remote control for easy and user-friendly operation. The units are equipped with a pressure pump optimized for internal circulation adjustable from performance level five to eight. To prevent moisture in the atmosphere from condensing at low temperatures, bath bridge and bath edge heating are integrated into the design. Proline Kryomats stand out for having the latest technologies and an excellent price-performance ratio.

<mark>ب</mark>

Application examples

Constant temperatures

- Notch bending test
- Drop test

Changing temperatures

- Determination of pour point
- Brookfield test of lubricants
- Test of slide bearings

Your advantages at a glance

+	The Proline Kryomats advantages	Your benefits
	 Removable Command remote control with graphic LCD Automatic adjustment of the control parameters via integrated software for adaptive control 	 Easy and intuitive operation. Quick setting changes Saves time-consuming calculation of control parameters
	 Offset control head Integrated bath edge and bath bridge heating Use of innovative cooling technology 	 Allows installation of optional supplementary pumps for external applications Avoids condensation and ice build-up High cooling capacity and low operating temperatures with very small footprint
	 Updated, adjustable pump nozzle 	 Optimum circulation and temperature distribution throughout the entire bath
	 Spacious baths with large bath openings Thread sleeves as standard on the edge of the bath 	 Accomodates various sample shapes and sizes with efficient flow Allow for the fixing of testing equipment without further conversion measures
	 Intelligent cooling fan control Optimised cooling airflow Internal release valve 	 Optimum heat discharge while reducing noise emission Bath drain at front of unit No protruding release valve

LAUDA Proline Kryomats

Proline Kryomats Air-cooled cooling thermostats

The air-cooled Proline Kryomats have a working temperature range from -50 and -90 up to 200 °C. The models are available with bath volumes of 30 and 40 liters. The Proline SmartCool system, with its energysaving digital cooling management, ensures that the cooling output is run in accordance with the application needs. That saves up to 75 percent of energy compared to standard cooling methods. Two different booster pumps are available as options (ex works) especially for external applications that require a considerable increase in volume flow/ discharge pressure.



Cooling curves Heat transfer liquid: Ethanol, bath closed Bath temperature 20 0 ① RP 3050 C 2 RP 4050 C -20 3 RP 3090 C ④ RP 4090 C -40 2 -60 1 4 -80 -100 0 10 20 30 40 50 60 70 80 90 100 min Cooling time

Temperature range -90...200 °C

*

Included accessories Bath cover · 4 closing plugs for pump connections · 2 connectors 13 mm

Additional accessories

Interface modules: analog, RS 232/485, contact, Profibus module

Options Booster pumps


		- 20-			
All technical data on page 96 a Other power supply variants on	•		1160 mm		
Technical features		RP 3050 C	RP 4050 C	RP 3090 C	RP 4090 C
Working temperature range*	°C	-50200	-50200	-90200	-90200
Temperature stability	±Κ	0.05	0.05	0.05	0.05
Heater power	kW	3.5	3.5	3.5	3.5
Cooling output at 20 °C	kW	5.0	5.0	3.0	3.0
Pump pressure max.	bar	0.5	0.5	0.5	0.5
Pump flow (pressure) max.	L/min	19	19	19	19
Bath volume	L	2331	3244	2331	3244
Bath opening/depth	mm	350x200/250	350x350/250	350x200/250	350x350/250
Cat. No. 400 V; 3/N/PE; 50 Hz		LUK 239	LUK 241	LUK 245	LUK 247

* Working temperature range is equal to ACC range

Proline Kryomats Water-cooled cooling thermostats

In the case of the water-cooled Proline Kryomats, the process heat is dissipated with the use of facility cooling water. This largely prevents unnecessary heating of the surrounding environment. As a result of this type of cooling, even higher cooling capacities are achieved than with the aircooled units. The electronic cooling water management minimizes water consumption. The booster pumps, available as options (ex works), are particularly recommended for external applications where increased volume flow or greater pressures are required.





Cooling curves Heat transfer liquid: Ethanol, bath closed

① RP 3050 CW 2 RP 4050 CW ③ RP 3090 CW ④ RP 4090 CW Applications Advantages Devices Accessories

Temperature range -90...200 °C

*

Included accessories

Bath cover · 4 closing plugs for pump connections · G ³/₄" locknut with 1/2" hose clip · 2 connectors 13 mm

Additional accessories

Tubing for cooling water · Interface modules: analog, RS 232/485, contact, Profibus module

Options Booster pumps

Cooling thermostat RP 4090 CW	STAR				
All technical data on pa Other power supply var		1160 mm			
Technical features		RP 3050 CW	RP 4050 CW	RP 3090 CW	RP 4090 CW
Working temperature range*	°C	-50200	-50200	-90200	-90200
Temperature stability	±Κ	0.05	0.05	0.05	0.05
Heater power	kW	3.5	3.5	3.5	3.5
Cooling output at 20 °C	kW	6.0	6.0	4.0	4.0
Pump pressure max.	bar	0.5	0.5	0.5	0.5
Pump flow (pressure) max.	L/min	19	19	19	19
Bath volume	L	2331	3244	2331	3244
Bath opening/depth	mm	350x200/250	350x350/250	350x200/250	350x350/250

LUK 242

LUK 246

LUK 248

LUK 240

* Working temperature range is equal to ACC range

Cat. No. 400 V; 3/N/PE; 50 Hz

Proline accessories

Shut down valve/Reverse flow protection Reverse flow protection when thermostating external systems, to avoid over-flow when pump stops, for retrofitting with LiBus. Temperature range -40...140 °C

Cat. No.	Description
LCZ 9673	Shut down valve reverse flow protection with LiBus
Suitable for	All Proline devices

Solenoid valve

Water-conscious cooling on heating thermostats for cooling water control. Controlled cooling operation for exothermal reactions or controlled cooling with programmer. Up to 155 $^{\circ}$ C bath temperature.

Cat. No.	Description	Temperature range
LCZ 9662	Solenoid valve with LiBus-connector	-10155 °C
Suitable for	All heating and clear-view thermostats	

Baskets

For notch bending test

Cat. No.	Suitable for
LCZ 0658	RP 870, RP 870 C, RP 890, RP 890 C
LCZ 0694	RP 1290, RP 1290 C

Constant level device

Necessary for the constant liquid level when thermostating open external baths. Connection set: for wall thickness of bath vessel between 0 to 30 mm, with opening for thermometers 4 mm or 1.9 mm \emptyset and clamping gland HX 077 and HX 078.

Cat. No.	Description	Suitable for
LCZ 0660	Level controller, mechanical	P 8 (C), RP 845 (C) RP 855 (C)*, RP 870 (C)*
LCZ 0679	Connection set for external inlet and outlet	LCZ 0660
* Not with option I	path cover including bath edge heating (LCZ 9670)	

Automatic filling device

For automatic replacement of liquid losses in thermostat bath, for example by evaporation. Also from vessels with max. 1 m suction height

Cat. No.	Description
LCZ 9661	Automatic filling device with LiBus

Controlled high-temperature chiller HTC with LiBus For controlled cooling of thermostats in the operating temperature range up to 300 °C without formation of vapors, to be connected to external water cooling source.

Cat. No.	Description	
LCZ 9663	Controlled high-temperature chiller HTC	



LCZ 9673



LCZ 9662



LCZ 0694









LCZ 9661

Proline Kryomats accessories

Interface modules

An RS 232/485 interface is integrated as a standard feature. The control head is equipped for two interface modules to be plugged into the rear of the unit.

Cat. No.	Description
LRZ 912	Analog module, 2 x In, 2 x Out, 0(4)20 mA or 010 V
LRZ 913	RS 232/485 interface, electrically isolated, 9-pin SUB-D socket
LRZ 914	Contact module NAMUR, 1 x In, 1 x Out, NE 28, 2 DIN socket
LRZ 915	Contact module SUB-D, 3 x In, 3 x Out, 15-pin SUB-D
LRZ 917	Profibus module, electrically isolated, 9-pin SUB-D socket

Suitable hoses/tubing for heat transfer liquids and cooling water Available upon request.

Booster pumps (only ex works)

For higher flow rates and pressure for external systems, connections M30 x 1.5 O

Cat. No.	Temperature range	Pressure max.	Pump flow max.
LWZ 080	-100150 °C	0.9 bar	90 L/min
LWZ 086	-40150 °C	3.2 bar	40 L/min
(O = outer thread)	1		

Baskets

For notch bending test

Cat. No.	Suitable for
LUZ 008	RP 3050 C, RP 3050 CW, RP 3090 C, RP 3090 CW
LUZ 009	RP 4050 C, RP 4050 CW, RP 4090 C, RP 4090 CW

Pour point determination

Bath cover accomodates up to 16 metal beakers

Cat. No.	Suitable for
UP 065	RP 4050 C, RP 4050 CW, RP 4090 C, RP 4090 CW



Order the detailed LAUDA accessories brochure and the heat transfer liquids brochure free of charge. These and additional product information can also be found at **www.lauda.de**





LRZ 912 LRZ 913 LRZ 914 LRZ 915 LRZ 917





LWZ 080



LUZ 008



UP 065