



Spray bottle with pump vaporizer



LaboPlast

Spray bottles with pump vaporizer and safety cap for the pump. Liquids can be finely vaporized at the push of a button in an environmentally compatible way. The transparent bottles make identification of contents easy. The protective cover prevents unintentional pressing of the pump and acts as a dust cover.

- Bottle made of HDPE, spray pump made of PP and stainless steel

Spray bottle with pump vaporizer

Capacity ml (oz.)	Dia. mm (in.)	Height with cap mm (in.)	Item no.
20 (0.68)	28 (1.10)	85 (3.35)	0309-3002 N3
50 (1.69)	35 (1.38)	110 (4.33)	0309-3005 N3
100 (3.38)	37 (1.46)	161 (6.34)	0309-3100 N3
250 (8.45)	50 (1.97)	222 (8.74)	0309-3025 N3