



Spray bottle



Pressureless spray bottle made of LDPE with long-life hand pump. The functioning of each individual sprayer is tested during manufacture.

Areas of application: Industry, laboratories, production, leisure, household. Suitable for agricultural pesticides and detergents.

- ▶ Media-contacting parts PP, PE and stainless steel
- ▶ Nozzle 0.6 mm (0.02 in.) diameter, nozzle steplessly adjustable
- ▶ Spray volume per stroke: 1.2 ml \pm 0.1 (0.04 oz. \pm 0.003)



Spray bottle

Capacity ml (oz.)	Height with cap mm (in.)	Neck dia. inside mm (in.)	Item no.
250 (8.45)	220 (8.66)	18 (0.71)	0309-1025 N4
500 (16.91)	240 (9.45)	18 (0.71)	0309-1050 N4
1000 (33.81)	295 (11.61)	17.5 (0.69)	0309-1100 N4



Spray bottle Turn 'n' Spray | Spray bottle with overhead valve



- ▶ Uninterrupted spraying, even over your head
- ▶ Parts in contact with media made of PP, PE, silicon and stainless steel
- ▶ Nozzle diameter 0.6 mm (0.02 in.), nozzle is continuously adjustable
- ▶ Amount sprayed per stroke: 1.2 ml ± 0.1 (0.04 oz. ± 0.003)

Spraying over your head is not possible with standard spray bottles. As soon as the suction tube is no longer immersed in liquid, you can't spray any more.

No so with our overhead spray bottle Turn'n'Spray. It features a special technology that allows you to spray from any position. That means you can spray practically and conveniently from any direction and reach hard-to-access places easily and without straining.

The high-quality, refillable quality spray bottle can be used again and again.

Spray bottle Turn 'n' Spray

Capacity ml (oz.)	Height with cap mm (in.)	Neck dia. inside mm (in.)	Item no.
250 (8.45)	220 (8.66)	18 (0.71)	0309-0002 N4
500 (16.91)	240 (9.45)	18 (0.71)	0309-0005 N4
1000 (33.81)	295 (11.61)	17.5 (0.69)	0309-0010 N4



Spray bottle with pump vaporizer



LaboPlast

Spray bottles with pump vaporizer and safety cap for the pump. Liquids can be finely vaporized at the push of a button in an environmentally compatible way. The transparent bottles make identification of contents easy. The protective cover prevents unintentional pressing of the pump and acts as a dust cover.

- ▶ Bottle made of HDPE, spray pump made of PP and stainless steel

Spray bottle with pump vaporizer

Capacity ml (oz.)	Dia. mm (in.)	Height with cap mm (in.)	Item no.
20 (0.68)	28 (1.10)	85 (3.35)	0309-3002 N3
50 (1.69)	35 (1.38)	110 (4.33)	0309-3005 N3
100 (3.38)	37 (1.46)	161 (6.34)	0309-3100 N3
250 (8.45)	50 (1.97)	222 (8.74)	0309-3025 N3