

Special waste containers in PP

5 litres capacity



SQUARE FORM.

Dimensions mm: 167x167x372.

A complete line of needle, sharp, pungent waste and hospital wastes, made of ecological plastic material (color-coated polypropylene copolymer complying with REACH, cadmium-free and other heavy metals).

They are resistant to fluids and solvents and are completely disposable in the incinerator without the release of harmful gases.

The entire line is designed to withstand perforation by needles and cannulae, the main cause of injuries between hospital staff.

Containers comply with applicable European and International standards. Depending on the type of container, they have NF302 Quality Certificate in accordance with NF X 30-500, conforming to UN-Certified ADR Standards and are Kitemark® certified according to the FDA-approved British Standard BS 7320 standard.

Autoclavable up to 134 ° C for 18 min.

Six measurements of square-shaped containers, with swivel cover that allows easy maneuverability for one-handed, temporary opening, and a simple, practical, snap-on closure.

Indication of maximum fill level.