

TERMOMETRI USO CHIMICA - SCALA OPALE - A LIQUIDO ECOLOGICO - Immersione parziale
GENERAL PURPOSE THERMOMETERS, ENCLOSED SCALE, LIQUID FILLING, PARTIAL IMMERSION

ART.	GRAD./RANGE °C	DIV.
100 101	-30 + 50	1
100 102	-10 + 50	1
100 103	-10 + 110	1
100 104	-10 + 150	1
100 105	-10 + 200	1
100 106	-10 + 250	1
100 110	-10 + 50	0,5
100 111	-10 + 100	0,5
100 112	-10 + 150	0,5
100 113	-10 + 200	0,5
100 114	-10 + 250	0,5



100 101

TERMOMETRI USO CHIMICA - A LIQUIDO ROSSO - A BASSO COSTO
 SCALA INCISA esternamente su vetro pieno fondo giallo - immersione totale - **Minimo 30 pezzi per tipo**
 GENERAL PURPOSE THERMOMOMETERS, RED LIQUID, LOW COST- TOT. IMMER.- SMALLEST 30 P.CES

ART.	GRAD./RANGE °C	DIV.
G11 582	-10 +110	1
G11 586	-10 +150	1



G11 582