

## CONTENITORE PER TRASPORTO VACCINI RCW12



<b>Volume lordo</b>	24 L   0,85 piedi cubi
<b>Capacità di stoccaggio del vaccino</b>	7 L
<b>Classe climatica (intervallo di temperatura ambiente)</b>	Zona calda (+ 43 ° C)
<b>Certificazioni</b>	CHI PQS E004 / 004
<b>Dimensioni esterne (AxLxP)</b>	499 x 550 x 475 mm   19,65 x 21,65 x 18,70 pollici
<b>Dimensioni interne (AxLxP)</b>	270 x 340 x 260 mm   10,63 x 13,39 x 10,24 pollici
<b>Peso netto</b>	11,7 kg   25,79 libbre
<b>Peso della Spedizione</b>	17 kg
<b>Materiali</b>	Polietilene
<b>Monitoraggio remoto</b>	Non remoto   tramite Data logger
<b>Vita a freddo a + 43 ° C</b>	114h54 min fino a 10 ° C
<b>Codice PQS dell'OMS</b>	E004 / 004
<b>Tipologia di prodotto</b>	Contenitore per il trasporto del vaccino

### External dimensions

Height (mm)	499
Width (mm)	550
Depth (mm)	475

### Internal dimensions

Height (mm)	270
Width (mm)	340
Depth (mm)	260
Storage height (mm)	209
Storage width (mm)	251
Storage depth (mm)	176

### Volumes

Gross volume (litres)	24
Storage volume (litres)	9.23

### Ship dimensions

Weight net (kg)	11.7
Height (mm) Cardboard	530
Width (mm) Cardboard	570
Depth (mm) Cardboard	490
Volume (m3) Cardboard	0.148
Weight (kg) Cardboard	17
Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

### Accessories

14 Ice Packs of 0.6L crystal thermometer
--

### Option

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

### Passive

### Vaccine transport box

### RCW12

### Technical details

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	90-105 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	23.3
Required Ice Packs	14 x 0.6L
Required Eutectics	-

### Performance

Hold over time at -10°C	-
Hold over time at +22°C	-
Hold over time at +32°C	-
Hold over time at 43°C	-
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	-
Cold life at +43°C	114.9h
Warm life at -20°C	40.9h
Cool life at +43°C	26.4h
Vaccine storage capacity (litres)	7
Climate Zone	-

### Certification

WHO test procedure	E4 / CB01-VP.1
WHO specifications	E4 / CB01.1
WHO test report	Intertek A65524c
PQS Code	E004/004
PIS Code	-
ADR/RID MARITIM IATA / ICAO	UN 4H2 / X 29.6 / S / "YYYY" B / 921-080154