

## CONTENITORE PER TRASPORTO VACCINI RCW8



<b>Volume lordo</b>	20 L   0,71 piedi cubi
<b>Capacità di stoccaggio del vaccino</b>	6 L
<b>Classe climatica (intervallo di temperatura ambiente)</b>	Zona calda (+ 43 ° C)
<b>Certificazioni</b>	CHI PQS E004 / 003
<b>Dimensioni esterne (AxLxP)</b>	437 x 588 x 288 mm   17,20 x 23,15 x 11,34 pollici
<b>Dimensioni interne (AxLxP)</b>	245 x 460 x 180 mm   5,71 x 18,11 x 7,09 pollici
<b>Peso netto</b>	6,8 kg   14,99 libbre
<b>Peso della Spedizione</b>	11 kg
<b>Materiali</b>	Polietilene
<b>Monitoraggio remoto</b>	Non remoto   tramite Data logger
<b>Vita a freddo a + 43 ° C</b>	57h54 min fino a 10 ° C
<b>Codice PQS dell'OMS</b>	E004 / 003
<b>Tipologia di prodotto</b>	Contenitore per il trasporto del vaccino

**External dimensions**

Height (mm)	437
Width (mm)	588
Depth(mm)	288

**Internal dimensions**

Height (mm)	245
Width (mm)	460
Depth (mm)	180
Storage height (mm)	202
Storage width (mm)	326
Storage depth (mm)	107

**Volumes**

Gross volume (litres)	20
Storage volume (litres)	7.05

**Ship dimensions**

Weight net (kg)	6.8
Height (mm) Cardboard	470
Width (mm) Cardboard	610
Depth (mm) Cardboard	310
Volume (m3) Cardboard	0.089
Weight (kg) Cardboard	11
Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

**Accessories**

10 Ice Packs of 0.6L 2 Ice Packs of 0.3L crystal thermometer
--

**Option**

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

<b>Passive</b>
<b>Vaccine transport box</b>
<b>RCW8</b>

**Technical details**

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	50-60 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	16.4
Required Ice Packs	10 x 0.6L and 2 x 0.3L
Required Eutetics	-

**Performance**

Hold over time at -10°C	-
Hold over time at +22°C	-
Hold over time at +32°C	-
Hold over time at 43°C	-
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	-
Cold life at +43°C	57.9h
Warm life at -20°C	21.6h
Cool life at +43°C	12h
Vaccine storage capacity (litres)	6
Climate Zone	-

**Certification**

WHO test procedure	E4 / CB01-VP.1
WHO specifications	E4 / CB01.1
WHO test report	Intertek A65524b
PQS Code	E004/003
PIS Code	-
ADR/RID MARITIM IATA / ICAO	UN 4H2 / Y19 / S / "YYYY" B / 921-110058