

Liofilchem®

CONTACT SLIDE

flexible contact slide for microbiological screening of surfaces

ISO 18593



Description	Packaging	Ref.
CONTACT SLIDE Chrom 1 (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides 120 slides	525292 # 53529 #
CONTACT SLIDE Chrom 2 (Chromatic™ E.coli-coliform / PCA+TTC + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides 120 slides	525392 53539
CONTACT SLIDE Chrom 3 (Chromatic™ E.coli-coliform / BAIRD PARKER AGAR + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides 120 slides	525462 # 53546 #
CONTACT SLIDE 1 (PCA + TTC + NEUTRALIZING / V.R.B.G. AGAR + NEUTRALIZING) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526
CONTACT SLIDE 2 (PCA + TTC + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527
CONTACT SLIDE 3 (PCA+TTC + NEUTRALIZING / VOGEL JOHNSON AGAR) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530
CONTACT SLIDE 4 (CETRIMIDE AGAR + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 # 53528 #
CONTACT SLIDE 5 (V.R.B.G. AGAR + NEUTRALIZING / BILE ESCULIN AGAR + NEUTRALIZING) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531
CONTACT SLIDE 8 (PCA + TTC + NEUTRALIZING / V.R.B.L. AGAR + NEUTRALIZING) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562
CONTACT SLIDE 9 (PCA + TTC + NEUTRALIZING / PCA + TTC + NEUTRALIZING) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532
CONTACT SLIDE 10 (ROSE BENGAL CAF AGAR + NEUTRALIZING / VOGEL JOHNSON AGAR) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #
CONTACT SLIDE TSA + NEUTRALIZING / TSA + NEUTRALIZING Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525472 # 53547 #
CONTACT SLIDE TSA + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 # 53548 #